

# ABOUT COMPANY

Get to know us.



Arrow Med Was Founded in 2017,  
To Manufacture and Distribute medical  
and Safety Technology Products .  
in Doing So, We Protect,Support and  
Save People`s lives,and Export all  
around The World.





**1ml**

**2ml**

**3ml**

**5ml**

**10ml**

**20ml**

**50ml**

**60ml**

**PRODUCTS SIZE**



# ArrowMed Normal syringe

With Exchangeable needle

in Compliance With Criteria Given in The ISO 7886-1 -

Normal Syringe With Exchangeable Needle

Capacity	Gauge size	Outer diameter (mm)	Length		Box/ Pcs	Syringe Tips
			(mm)	(inch)		
1 ml	25G	0.5mm	25mm	1"	300 Pcs	Luer Slip Syringe
2 ml	23G	0.6mm	25mm	1"	100/300 Pcs	
3 ml	23G	0.6mm	40mm	1 1/2"	100/300 Pcs	Luer lock syringe
	23G	0.6mm	32mm	1 1/4"		
	23G	0.6mm	25mm	1"		
	24G	0.55mm	32mm	1 1/4"		
	24G	0.55mm	25mm	1"		
	25G	0.5mm	25mm	1"		
5 ml	22G	0.7mm	32mm	1 1/4"	100/200 Pcs	Luer lock syringe
	23G	0.6mm	32mm	1 1/4"		
	23G	0.6mm	25mm	1"		
10 ml	21G	0.8 mm	25mm	1"	100/150 Pcs	Luer lock syringe
	21G	0.8 mm	40mm	1 1/2"		
20 ml	18G	1.2 mm	40mm	1 1/2"	100 Pcs	Luer lock syringe
	21G	0.8 mm	40mm	1 1/2"		
50 ml	18G	1.2 mm	40mm	1 1/2"	50 Pcs	Luer Slip Syringe
	21G	0.8 mm	40mm	1 1/2"		
60 ml	18G	1.2 mm	40mm	1 1/2"	50 Pcs	Luer Slip Syringe
	21G	0.8 mm	40mm	1 1/2"		



# ArrowMed Auto Disable Syringe

With Exchangeable needle

in Compliance With Criteria Given in The **ISO 7886-4** -  
Prevention From Sharp injuries in The Hospital and Healthcare sector



open The cover as Shown  
in This Pictures



Push The plunger Thoroughly To  
The End Of Complete injection



Push Plunger To Ensure That The Self-Closure  
Mechanism Working Well

## Auto Disable Syringe With Exchangeable Needle

Capacity	Gauge size	Outer diameter (mm)	Length		Box/ Pcs
			(mm)	(inch)	
1 ml	25G	0.5mm	25mm	1"	300 Pcs
2ml	23G	0.6mm	25mm	1"	100/250 Pcs
3ml	23G	0.6mm	40mm	1 1/2"	100 Pcs
	23G	0.6mm	32mm	1 1/4"	
	23G	0.6mm	25mm	1"	
	24G	0.55mm	32mm	1 1/4"	
	24G	0.55mm	25mm	1"	
	25G	0.5mm	25mm	1"	
5 ml	25G	0.5mm	16mm	5/8"	100/200 Pcs
	22G	0.7mm	32mm	1 1/4"	
	23G	0.6mm	32mm	1 1/4"	
	23G	0.6mm	25mm	1"	
10 ml	21G	0.8 mm	25mm	1"	100 Pcs
	21G	0.8 mm	40mm	1 1/2"	



Luer Slip Syringe




Luer lock syringe

# ArrowMed INSULIN Syringe

in Compliance With Criteria Given in The **ISO 8537**

Capacity	Gauge size	Outer diameter (mm)	Length		Box/ Pcs	
			(mm)	(inch)		
1 ml(100 U)	<b>28 G</b>	0.37mm	13mm	1/2"	100/300 Pcs	

# ArrowMed ORAL Syringe (Non-Sterile)

Capacity	Box/ Pcs	
5 ml	100 PCS	
10 ml	80 pcs	

# ArrowMed Catheter Syringe (Non-Sterile)

Capacity	Box/ Pcs	
50 ml	50 pcs	

# GET IN TOUCH

+2 01029998381  
+2 01096993522



+2 088-496-4893



<https://arrow-med.com>



[info@arrow-med.com](mailto:info@arrow-med.com)



ARROW FOR MEDICAL NECESSITIES

Arrow MED



ARB EL AWAMER INDUSTRIAL ZONE  
ABNUB, ASSUIT, EGYPT



[\Arrowmed2015](#)